POSTER PRESENTATION



Open Access

Dynamic cyclin profiles as a tool to segregate the cell cycle

David Garcia Munzer¹, Margaritis Kostoglou², Michalis C Georgiadis³, Efstratios N Pistikopoulos¹, Athanasios Mantalaris^{1*}

From 23rd European Society for Animal Cell Technology (ESACT) Meeting: Better Cells for Better Health Lille, France. 23-26 June 2013

Background and novelty

Mammalian cells growth, productivity and cell death are highly regulated and coordinated processes. The cell cycle is at the centre of cellular control and should play a key role in determining optimization strategies towards improving productivity [1]. Specifically, cell productivity is cell cycle, cell-line and promoter dependant [2]. The cyclins are key regulators that activate their partner cyclin-dependent kinases (CDKs) and target specific proteins driving the cell cycle. To our knowledge, there is no information on cyclin phase-dependent expression profiles of industrial relevant mammalian cell lines. We use the cyclin profiles as a tool to identify and quantify the landmarks of the cell cycle and implement a modelling approach to describe the bioprocess. Hereby, we introduce two possible experimental approaches to obtain such dynamic cyclin profiles.

Experimental approach

Cyclin expression (cyclin E - G_1 class and cyclin B - G_2 class) was studied in GS-NS0 batch cultures by flow cytometry. Two set of experiments were performed: a) culture of cells under perturbed (cell arrest) and unperturbed growth (control run) and b) culture of cells for DNA labelling to perform a proliferation assay as well as a non-exposed cells (control run). The static profiles were obtained by direct cyclin staining and the dynamic profiles were reconstructed by either a) tracking a partially synchronized population or b) combining the timings from proliferation assays with the static profiles.

Result discussion

Both cyclins showed a clear cell cycle phase-specific pattern (cyclin E was 10% higher at G_1 and cyclin B was 40% higher at G_2). These results were consistent among all the different culture conditions and were inferred from the static cyclin profiles. After the arrest release the dynamic cyclin profiles can be directly reconstructed by plotting the relevant cyclin content from the partially synchronized moving population traversing the cycle. An advantage of this approach is a clear view of the cyclin accumulation and transition threshold levels. However, this approach requires testing using different arrest agents, exposure levels and timings, which could have an effect on the cell behaviour.

A second approach included an indirect dynamic cyclin profile reconstruction by combining the acquired proliferation times for different cell cycle phases (e.g. G_1/G_0 , G_2/M) with the static cyclin profiles. If the static cyclin profiles are considered as the most representative cyclin values (and near to the transition threshold level), it is possible to reconstruct the dynamic profile by linking the threshold values with the cycling times (from the proliferation assay). The advantage of such approach is the ability to formulate different dynamic cyclin profiles such as constant functions, piece-wise linear functions or more elaborated profiles. However, implementation of such an approach requires the tuning of the proliferation assay and the frequency of sampling since it will affect the quality of the assay.

The two approaches showed comparable results both for the static cyclin profiles (also when compared to the control runs) and the dynamic cyclin profiles.

Conclusions

The different approaches for deriving the dynamic cyclin profiles provide a versatile experimental toolbox for cell



© 2013 Munzer et al.; licensee BioMed Central Ltd. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/2.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited. The Creative Commons Public Domain Dedication waiver (http:// creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.

^{*} Correspondence: amantalaris@imperial.ac.uk

¹Biological Systems Engineering Laboratory, Centre for Process Systems Engineering, Department of Chemical Engineering, Imperial College London, London, SW7 2AZ, UK

Full list of author information is available at the end of the article

cycle characterization. Cyclins can be used as cell cycle distributed variables and be experimentally validated (quantitatively), avoiding the use of weakly supported variables (e.g. age or volume). The observed patterns and timings provide a blueprint of the cell line's cell cycle, which can be used for cell cycle modelling. The development of these models will aid the systematic study of the cell culture system, the improvement of productivity and product quality.

Acknowledgements

The authors are thankful for the financial support from the MULTIMOD Training Network, European Commission, FP7/2007-2013, under the grant agreement No 238013 and to Lonza for generously supplying the GS-NS0 cell line.

Authors' details

¹Biological Systems Engineering Laboratory, Centre for Process Systems Engineering, Department of Chemical Engineering, Imperial College London, London, SW7 2AZ, UK. ²Department of Chemistry, Aristotle University of Thessaloniki, Thessaloniki, 54124 Greece. ³Department of Chemical Engineering, Aristotle University of Thessaloniki, Thessaloniki, 54124 Greece.

Published: 4 December 2013

References

- 1. Dutton RL, Scharer JM, Moo-Young M: Descriptive parameter evaluation in mammalian cell culture. *Cytotechnol* 1998, **26**:139-152.
- Alrubeai M, Emery AN: Mechanisms and Kinetics of Monoclonal-Antibody Synthesis and Secretion in Synchronous and Asynchronous Hybridoma Cell-Cultures. J Biotechnol 1990, 16:67-86.

doi:10.1186/1753-6561-7-S6-P523

Cite this article as: Munzer *et al.*: Dynamic cyclin profiles as a tool to segregate the cell cycle. *BMC Proceedings* 2013 7(Suppl 6):P23.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

BioMed Central

Submit your manuscript at www.biomedcentral.com/submit