

POSTER PRESENTATION

Open Access

A fast solution to NGS library preparation with low nanogram DNA input

Pingfang Liu*, Gregory JS Lohman, Eric Cantor, Bradley W Langhorst, Erbay Yigit, Lynne M Apone, Daniela B Munafo, Christine Sumner, Fiona J Stewart, Thomas C Evans Jr, Nicole M Nichols, Eileen T Dimalanta, Theodore B Davis

From Beyond the Genome 2012 Boston, MA, USA. 27-29 September 2012

Next-generation sequencing (NGS) has significantly impacted human genetics, enabling a comprehensive characterization of human genome as well as better understanding of many genomic abnormalities. By delivering massive DNA sequences at unprecedented speed and cost, NGS promises to make personalized medicine a reality in the foreseeable future. To date, library construction with clinical samples has been a challenge, primarily due to the limited quantities of sample DNA available. To overcome this challenge, we have developed a fast library preparation method using novel NEBNext reagents and adaptors, including a new DNA polymerase that has been optimized to minimize GC bias. This method enables library construction from an amount of DNA as low as 5 ng, and can be used for both intact and fragmented DNA. Moreover, the workflow is compatible with multiple NGS platforms.

Published: 1 October 2012

doi:10.1186/1753-6561-6-S6-P26

Cite this article as: Liu et al.: A fast solution to NGS library preparation with low nanogram DNA input. BMC Proceedings 2012 6(Suppl 6):P26.

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at www.biomedcentral.com/submit





